


Machine Learning based Epileptic Seizure Prediction

Aditya Sriram S
Sathiya Murthi S
Duke Daffin T


Under the guidance of
Dr. D Sangeetha - Assistant Professor, MIT
Dr. S Umamaheswari - Assistant Professor, MIT

Block Diagram

Training



Bi-LSTM Model Architecture



Model Summary

Total params: 2,690,475
Trainable params: 2,690,347
Non-trainable params: 128

Layer (type)	Output Shape	Param #
=====		
conv2d_3 (Conv2D)	(None, 1278, 22, 32)	224

max_pooling2d_3 (MaxPooling2)	(None, 639, 11, 32)	0

batch_normalization_2 (Batch Normalization)	(None, 639, 11, 32)	128

conv2d_4 (Conv2D)	(None, 637, 10, 32)	6176

max_pooling2d_4 (MaxPooling2)	(None, 318, 5, 32)	0

batch_normalization_3 (Batch Normalization)	(None, 318, 5, 32)	128

conv2d_5 (Conv2D)	(None, 316, 4, 32)	6176

max_pooling2d_5 (MaxPooling2)	(None, 158, 2, 32)	0

flatten_1 (Flatten)	(None, 10112)	0

repeat_vector_1 (RepeatVector)	(None, 1, 10112)	0

bidirectional_2 (Bidirectional)	(None, 1, 64)	2597120

bidirectional_3 (Bidirectional)	(None, 64)	24832

dense_5 (Dense)	(None, 300)	19500

dense_6 (Dense)	(None, 100)	30100


dense_7 (Dense)	(None, 50)	5050

dense_8 (Dense)	(None, 20)	1020

dense_9 (Dense)	(None, 1)	21
=====		

Block Diagram

Testing



Outcomes

- Machine Learning has been employed to classify raw EEG signal data into Pre-ictal and Inter-ictal periods.
- Classifier model is trained using CHB-MIT EEG Dataset.
- Classification is carried out on 5-minute intervals of 23-channel EEG signal collected using Brain-Computer Interface.
- The model makes use of Convolutional layers to find spatial relationships and Bi-directional LSTM layers to find temporal relationships, so as to classify pre-ictal and inter-ictal periods.
- When this model is fed real-time data, classification can be done in real-time so as to identify potential pre-ictal periods, allowing pre-emptive indication of an impending Epileptic seizure.
- This alarm would allow for medication to be taken ahead of an Epileptic Seizure to reduce the effects of the seizure or even avoid it completely.